




Neuropixels workshop



Paris Neuro course 2019
Optical and electrical recordings in neuroscience

Neuropixels 1.0



IMEC - SEM image of a probe tip and electrodes

IMEC – large scale manufacture



Resources:

Get probes and more

<https://www.neuropixels.org/>

Wiki

<https://github.com/cortex-lab/neuropixels/wiki>

Slack channel

<https://neuropixelsgroup.slack.com>

UCL course slides


<https://www.ucl.ac.uk/neuropixels/courses/2019-course>

How to videos:


<https://billkarsh.github.io/SpikeGLX/#how-to-videos>

Recording with neuropixels 1.0 probes


Xilinx




1 probe per system;
ethernet connection




PXI



4 probes per module; up to 16 per PXI



Neuropixels probes




*Standard PXIe hardware
+
Neuropixels module*

Image from: <http://www.open-ephys.org/neuropixels>

Planning a recording – stimuli

<https://www.youtube.com/watch?v=8VdFf3egwfg>




J. Physiol. (1959) 148, 574–591

RECEPTIVE FIELDS OF SINGLE NEURONES IN THE CAT'S STRIATE CORTEX

BY D. H. HUBEL* AND T. N. WIESEL*

*From the Wilmer Institute, The Johns Hopkins Hospital and
University, Baltimore, Maryland, U.S.A.*

(Received 22 April 1959)





Planning a recording – simultaneous V1 and dLGN

Using the Allen Common Coordinate Framework

Andy Peters <https://github.com/cortex-lab/allenCCF>



`allen_ccf_npx`

Probe insertion: -3187 AP, -2400 ML, Angle: 78 deg horizontal



Neuropixels 1.0 – probe holders



3D printed holder



Nick Steinmetz

https://github.com/cortex-lab/neuropixels/wiki/Probe_handling

Dovetail holder



HHMI Janelia

<https://flintbox.com/public/project/52167/>

Spike sorting - installation

Hardware:

CUDA capable graphics card (4Gb memory)
> 16Gb memory ($\frac{1}{4}$ of the max recording file recommended)
Fast disk - SSD

Software:

CUDA toolkit (depends on matlab version)
Visual Studio Community (on Windows)
Matlab (>R2016)
Toolboxes (check documentation):
 Signal Processing
 Parallel computing
 Statistics and Machine learning

Matlab version	CUDA toolkit
R2019a	10.0
R2018b	9.1
R2018a	9.0
R2017b	8.0
R2017a	8.0
R2016	7.5

Kilosort2 - <https://github.com/MouseLand/Kilosort2>

Download: `git clone https://github.com/MouseLand/Kilosort2.git`

Matlab install: `cd CUDA % in the kilosort2 folder`
`mexGPUall`

JRClust - <https://github.com/JaneliaSciComp/JRCLUST>

Download: `git clone https://github.com/JaneliaSciComp/JRCLUST.git`

Matlab install: `addpath('/path/to/JRCLUST')`
`jrclust.CUDA.compileCUDA();`

Spike sorting – running

Kilosort2 - <https://github.com/MouseLand/Kilosort2>

With the GUI: `kilosort2`

<https://github.com/MouseLand/Kilosort2/wiki/1.-The-GUI>

Through the command line:

- 1) **Create a channelmap:** look at `\configFiles\createChannelMap.m`
- 2) **Copy** the `master_kilosort.m` and `configFiles\StandardConfig_MOVEME.m` to another folder (e.g. data folder).
- 3) **StandardConfig.m:** **Set the paths** `ops.fbinary`, `ops.fproc`, `ops.chanMap` and `ops.root`, the **sampling frequency** – `ops.fs`, and the **number of channels** in the file – `ops.NchanTOT`.
- 4) `master_file.m`: Set the paths to the config file and the repositories.

Run the master_file.

JRClust - <https://jrclust.readthedocs.io/en/latest/usage/index.html>

Create a config file: `jrc bootstrap`

Display traces: `jrc traces <configfile.prm>`

Do the sorting: `jrc detect-sort <configfile.prm>`

Manual sorting GUI: `jrc manual <configfile.prm>`

Spike sorting – manual curation kilosort2


Installation:

Install dependencies using conda: <https://github.com/kwikteam/phy>

`pip install git+https://github.com/kwikteam/phy git+https://github.com/kwikteam/phy-contrib --upgrade`

Documentation:

<https://phy-contrib.readthedocs.io/en/latest/template-gui/> - Important! Has also decision criteria.




Spike sorting – manual curation JRClust

jrc manual <configfile.prm>

Documentation:

<https://jrclust.readthedocs.io/en/latest/pipeline/curate/index.html#the-waveform-view-figwav>



Spike sorting – manual curation

Documentation:

<https://phy-contrib.readthedocs.io/en/latest/template-gui/#a-typical-approach-to-manual-clustering>

Small waveforms?

Yes: label as MUA; move to the next.

No: compare with most similar waveforms.

Are the clusters in different channels?

Yes: Move to the next comparison.

No: inspect cross-correlograms.

Are the autocorrelograms different and there is no refractory period on the cross-correlogram?

Yes: Not the same neuron

No: Check amplitude stability and overlap in feature space

Units drift continuously in time and don't overlap?

Yes: Probably the same unit, merge.


No: Move to the next comparison.

Clear refractory period and all reasonable comparisons made?

Yes: Classify as “good”

No: Classify as “multi-unit activity” (MUA)

Neuropixels 3A – xilinx and basestation



External trigger (5V)
or start signal

